

Angelantoni Life Science

Refrigerators
Ekozero Series



EKOZERO

 **HeraScientific**
L I F E S C I E N C E

- Maximum uniformity of temperature
- Automatic defrost
- Forced air circulation
- Easy to clean
- Door with triple anti-moisture glass or solid door
- Alarm for highest/lowest temperature and door open
- Two capacities available: 700 lt and 1500 lt

ALS

→ Temperature range 0/+15°C

Technical features

Options

- Temperature chart recorder
- 23-50 mm port hole
- PT100 probe



STRUCTURE / Grey pre-painted metal sheet frame. Adjustable feet for proper positioning. Internal chamber in stainless steel AISI 304, with rounded corners to allow an easy cleaning.

DOOR / Reversible service door made with same features of the frame or in anti-moisture triple glass, equipped with a magnetic lock. Key lock to preclude access to unauthorized personnel. Mechanical door lock at 90°, to facilitate the loading and unloading of samples.

SHELVES / No. 3, in stainless steel, adjustable in height.

INNER LAMP / Automatic or manual courtesy light.

INSULATION / The 75 mm thickness offers an excellent insulation and a significant energy saving.

COOLING SYSTEM / Fully sealed, placed on the top side of the cabinet, featuring a low noise hermetic compressor and air condenser. Forced ventilation. Automatic or manual (in case of emergency) defrost.

REGULATION AND CONTROL SYSTEM:

- Luminous digital display
- Chamber On/Off switch
- ON/OFF manual defrost switch
- Working temperature settings
- ON/OFF internal light switch
- Acoustic and visual high/low temperature alarm
- Acoustic and visual open door alarm
- Highest/lowest temperature reached
- Predisposition for remote alarm contact NA 230 V

Model	EKOZERO			
	700/1 TN	700/1 TN - GL	1500/1 TN	1500/1 TN - GL
Capacity (l)	600		1400	
Temperature range (°C)	0 / +15			
Supply voltage (V)	230V(+6% -10%) / 50 (60) Hz / 1+ T(G)			
Maximum absorbed current (A)	3,3		3,3	3,5
Refrigerant	R404A			
External dimensions (mm) (WxDxH)	740x815x1998		1480x815x1998	
Internal dimensions (mm) (WxDxH)	590x675x1500		1330x675x1500	
Overall weight with standard packaging (Kg)	140	150	230	240
Noise level (dba)	49			

→ Illustration, description and technical data may be changed without notice.
 (*) 230 W/60 Hz models are also available



Herascientific Life Science S.L.

C/ Mistral Nº2, 139 - 28223 Pozuelo de Alarcón (Madrid)

Tel: + 34 91 679 99 59 - Mov: + 34 639 71 89 85

E-mail: ejecutivo@herascientific.com