



SICUR® 60

Technical data

CABINET FOR THE STORAGE OF DIFFERENT INFLAMMABLE LIQUIDS

- Completely realized in pickled steel 10/10 painted with dusts RAL 1021 colour
- Rounded off edges and invisible screws
- Air vents for natural air-recycle
- Metal lock with code pole and plastic/metal lock with for easy handle; supplied with 2 keys
- Reinforced steel hinges 15/10
- With 3 tray shaped shelves adjustable in height
- 1 bottom basin to collect eventual spread of liquids
- Earth ground wire
- Adjustable legs
- The cabinet is disassembled and supplied with an easy assembly kit.
- Warning stickers showing “Inflammable Products” and “No Smoking”
- Insurance RC products for all cabinets’ life.
- External dimensions (WxDxH) : mm 575 x 500 x 1850 h
- I Single packaging (pallet): mm 280 x 1860 x 600 h
- Weight: kg 57 (kg 67 with packaging)